OSYPKA iSense®

ELECTRICAL CARDIOMETRY™ SENSORS


Window to the Circulation®

OSYPKA iSense® sensors are intended for use with the Electrical Cardiometry™ (EC™) monitor AESCULON® and cardiometers ICON® and ICONCore®. OSYPKA iSense® Sensors are 100% disposable and are available in various sizes, for neonate, pediatric, and adults. All sensors utilize a highly conductive AgCl wet gel and strong hypoallergenic adhesion, which provides instant electrical contact for high quality signal conductivity. Four iSense® sensors are packed together in one pouch for the convenience of the user and more importantly keeping the gel fresh for best conductivity and signal quality, specially in difficult patients.

OSYPKA iSense® sensor box

Patient cable

OSYPKA iSense® Specification:

Sensor	Patient Group	Type/ Size	Article No.	Packaging Unit
	Premature and Pre-mature	Prewired sensor 22 mm x 22 mm (rectangular)	XN.TSP	20 pouches
	Adults and Children	Sensor with clip 30 mm (round)	XN.TSN	25 pouches

Products and options may not be available in every region or country. Please contact mail@osypkamed.com for availability in your region and country and for further information.

Osypka Medical GmbH

Albert-Einstein-Strasse 3
12489 Berlin, Germany
Phone: +49 (30) 6392 8300
Fax: +49 (30) 6392 8301
E-Mail: mail@osypkamed.com
www.osypkamed.com

United States of America: Osypka Medical, Inc.

dba. Cardiotronic
La Jolla, CA 92037, USA
Phone: +1 (858) 454 0600
Fax: +1 (858) 454 0640
E-Mail: mail@cardiotronic.net
www.cardiotronic.net

U.S. Patent Nr. 6,511,438. Other patents pending.

iSense, AESCULON, ICON, ICONCore, Cardiotronic, Electrical Cardiometry, Window to the Circulation and logos are trademarks of Osypka Medical.

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Vascular System

SVR / SVRI	Systemic Vascular Resistance / SVR-Index based on input of MAP and CVP
------------	--

Contractility

ICON™	Index of Contractility
VIC™	Variation of Index of Contractility
STR	Systolic Time Ratio
PEP	Pre-Ejection Period
LVT	Lvet-Ventricular Ejection Time
CPI	Cardiac Performance Index

Fluid Status

TFC	Thoracic Fluid Content
SVV	Stroke Volume Variation
FTC	Corrected Flow Time

Oxygen Status

DO2 / DO2I	Oxygen Delivery/DO2-Index based on input of Hemoglobin and SpO2
CaO2	Arterial Oxygen Content

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測量方法	Electrical Cardiometry (EC)
測量電流	≤ 2.0 mA RMS / 50 kHz
ECG	15...250 bpm (±5%/5 bpm)
電源輸入	100...240 VAC, 47...63 Hz
消耗功率	max. 15 VA
內置電池	Lithium Ion, cap. > 2 hours
顯示器	3.5" color TFT
外殼尺寸	205 mm x 110 mm x 38 mm
重量	750g
儀器類別	Class IIa · Class II equipment (Type BF)

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- Rajput R, Das S, Chauhan S, Bisoi A, Vasdev S. Comparison of Cardiac Output Measurement by Noninvasive Method with Electrical Cardiometry and Invasive Method with Thermodilution Technique in Patients Undergoing Coronary Artery Bypass Grafting. *World Journal of Cardiovascular Surgery*. 2014.
- Malik V, Subramanian A, Chauhan S, Hote M. Correlation of Electric Cardio-metry and Continuous Thermodilution Cardiac Output Monitoring Systems. *World Journal of Cardiovascular Surgery*. 2014.
- Peev M et al. Real-time sample entropy predicts life-saving interventions after the Boston Marathon bombing. *Journal of Critical Care*. 2013.
- Mejaddam A, King D, et al. Real-time heart rate entropy predicts the need for lifesaving interventions in trauma activation patients. *J Trauma Acute Care Surg*. 2013.
- Flinck et al. Cardiac output measured by electrical velocimetry in the CT suite correlates with coronary artery enhancement: a feasibility study. *Acta Radiol*. 2010.
- Zoremba et al. Comparison of electrical velocimetry and thermodilution techniques for the measurement of cardiac output. *Acta Anaesthesiol Scandinavia*. 2007.
- Schmidt et al. Comparison of electrical velocimetry and tranoesophageal Doppler echocardiography for measuring stroke volume and cardiac output. *British Journal of Anaesthesia*. 2005.

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- Lien R, Hsu KH, Chu J, Chang Y. Hemodynamic alterations recorded by electrical cardiometry during ligation of ductus arteriosus in preterm infants. *European Journal of Pediatrics*. 2014. 1-8.
- Coté CJ, et al. Continuous noninvasive cardiac output in children: is this the next generation of operating room monitors? Initial experience in 402 pediatric patients. *Paediatr Anaesth*. 2014.
- Grollmuss O et al. Non-invasive cardiac output measurement in low and very low birth weight infants: a method comparison. *Front Pediatr*. 2014.
- Noonan P, Viswanathan S, Chambers A, Stumper O. Non-invasive cardiac output monitoring during catheter interventions in patients with cavopulmonary circulations. *Cardiol Young*. 2014.
- Grollmuss O, Demontoux S, Capderou A, Serraf A, Belli E. Electrical velocimetry as a tool for measuring cardiac output in small infants after heart surgery. *Intensive Care Med*. 2012.
- Rauch R, Wlisch E, Lansdell N, et al. Non-invasive measurement of cardiac output in obese children and adolescents: comparison of electrical cardiometry and transthoracic Doppler echocardiography. *J Clin Monit Comput*. 2012.
- Noori S, Drabu B, Soleymani S, Seri I. Continuous Non-invasive cardiac output measurements in the neonate by electrical velocimetry: a comparison with echocardiography. *Arch Dis Child Fetal Neonatol* Ed 2012.
- Norozi K et al. Electrical velocimetry for measuring cardiac output in children with congenital heart disease. *Br J Anaesth*. 2007.
- Osthaus et al. Comparison of electrical velocimetry and transpulmonary thermodilution for measuring cardiac output in piglets. *Pediatric Anesthesia*. 2007.

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Copyright © 2024 Osypka Medical